RFC29400

**Indicazioni**

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RENOL products are safety, immobilization, restraint products

Indications

1. RENOL products are used for medical purposes to secure an agitated patient in the chair, to avoid an accident, to allow the patient to follow his medical treatment.

2. RENOL products are only used under medical prescription.

3. The products must be used by trained staff following the RENOL's instructions for safe use and according to the establishment's safety rules.

Contraindications

4. Ensure that RENOL products do not come into contact with injured skin, do not obstruct breathing nor impede blood pressure.

Precautions of use

5. To avoid any accident and in order to ensure the person's comfort and safety, it is essential to read the instruction leaflets carefully and understand them before operating belts. All users follow the instructions and procedures, patients' safety depends on using belts correctly.

6. Complications, accidents related to the restraint in case of misuse, during unsuitable, unsupervised use: bedsores, pneumonia or deep vein thrombosis, psychological stress (agitation, delirium) and negative impact on cognitive abilities, consequence fatal (following severe trauma or suffocation).

7. Ensure that a RENOL product is not damaged before use, all products must be clean and in good condition.

8. We recommend monitoring the immobilised patient frequently; restless and combative patients must be under constant observation.

9. Always use the correct size: too small or too large may compromise safety or comfort. Measure the patient's waist size, select the appropriate RENOL product.

10. Remove all sharp objects from the patient (jewellery, watches, glasses, etc.).

11. Only use products RENOL on suitable support. Never use RENOL products as a safety belt in vehicles.

12. When supporting or immobilising the patient in a chair, don't forget the patient's ability to turn over.

13. Connect RENOL products to a solid part of the chair. The straps must be secured; they cannot slide or change position; the fastening point of the strap must not damage the product (sharp edge, etc.).

14. Ensure that the RENOL key, part number: RFK00100 (RL) and RFK01111 (M) is always available. Use only RENOL locks.

15. In an emergency, if the key is not available or problem with lock system, cut the strap close to the frame.

16. Before use, check operation of the magnet (M) lock system, check correct operation of the RFK01202 Button with the key RFK01111 and pin RFK01101. We recommend replacing the Button RFK01202 after 2 years. The lifetime may be shorter or longer, depending on use.

17. Before use, check operation of the mechanical lock system (RL), check correct operation of the RFK00400 Button with the key RFK00100 and long pin RFK00300. We recommend replacing the Button RFK00400 every 3 years. The lifetime may be shorter or longer, depending on use.

18. The magnetic key (M) must be placed 10 cm minimum from a pacemaker or any other medical system that could be disturbed by a magnet.

19. The magnetic lock system could be activated with a magnet other than that delivered by RENOL.

Maintenance and materials

20. Follow the washing instructions marked on the products.

21. Close Velcro fastenings before washing and drying to prevent them wearing.

22. Detergents: avoid using chlorine.

23. Temperature: maximum temperature 95°C.

24. Drying: delicate cycle at low temperature, about 5% shrinkage possible.

25. Remove all locking systems before washing and drying.

26. Never immerse the magnetic (M) and mechanical closures (RL) in a liquid.

27. Material: 100% Terylene (polyester).

28. The raw materials are treated with fire retardant and anti bacterial.

29. 2 years warranty, this warranty does not cover damage resulting from normal wear and tear.

All our products are CE-approved - Class I medical devices.

Any serious incident occurring in connection with this device must be notified to the manufacturer and to the competent authority of the country where the serious incident took place.

For more information, contact your support specialist, your distributor or www.renolcare.com