N2

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Instructions for use ENGLISH



N2


REF	Désignation / Omschrijving Bezeichnung / Designation	Mesures / Maten / Maßnahmen / Measurements	REF	Désignation / Omschrijving Bezeichnung / Designation
RFA17200	STANDARD ANKLE HOLDERS M	18-28cm	RFK01111	MAGNET KEY
RFA17210	STANDARD ANKLE HOLDERS S	16-24cm	RFK01222	MAGNET BUTTON + PIN
RFA17220	STANDARD ANKLE HOLDERS L	23-32cm	RFK01201	MAGNET PIN
RFA17550	SECURITY ANKLE HOLDERS M	18-28cm	RFK01202	MAGNET BUTTON
RFL15000	LIMB HOLDERS	14-28cm	RFK00100	STEEL KEY
RFW16200	STANDARD WRIST HOLDERS M	14-22cm	RFK00300	STEEL PIN
RFW16210	STANDARD WRIST HOLDERS S	10-18cm	RFK00400	STEEL BUTTON
RFW16220	STANDARD WRIST HOLDERS L	18-26cm	RFK00433	STEEL BUTTON + PIN
RFW16550	SECURITY WRIST HOLDERS M	14-22cm		
			REF + M =	BELT + MAGNET LOCKS
			REF + RL =	BELT + STEEL LOCKS

REF + M =

REF + RL =

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RENOL products are safety, immobilization, restraint products

Indications

1. RENOL products are used for medical purposes to secure an agitated patient in bed, on a stretcher, to avoid an accident, to allow the patient to follow his medical treatment.
2. RENOL products are only used under medical prescription.
3. The products must be used by trained staff following the RENOL's instructions for safe use and according to the establishment's safety rules.

Contraindications

4. Ensure that RENOL products do not come into contact with injured skin, do not obstruct breathing nor impede blood pressure.

Precautions of use

5. To avoid any accident and in order to ensure the person's comfort and safety, it is essential to read the instruction leaflets carefully and understand them before operating belts. All users follow the instructions and procedures, patients' safety depends on using belts correctly.
6. Complications, accidents related to the restraint in case of misuse, during unsuitable, unsupervised use: bedsores, pneumonia or deep vein thrombosis, psychological stress (agitation, delirium) and negative impact on cognitive abilities, consequence fatal (following severe trauma or suffocation).
7. Ensure that a RENOL product is not damaged before use, all products must be clean and in good condition.
8. We recommend monitoring the immobilised patient frequently; restless and combative patients must be under constant observation.
9. Always use the correct size: too small or too large may compromise safety or comfort. Measure the patient's waist size, wrist and ankle, select the appropriate RENOL product.
10. Remove all sharp objects from the patient (jewellery, watches, glasses, etc.).
11. Only use products RENOL on suitable support. Never use RENOL products as a safety belt in vehicles.
12. When supporting or immobilising the patient in a bed, don't forget the patient's ability to turn over.
13. Connect RENOL products to a solid part of the bed, stretcher. The straps must be secured; they cannot slide or change position; the fastening point of the strap on the bed must not damage the product (sharp edge, etc.).
14. The straps (RFW16550-RFA17550) are placed on a part of the bed that can move to match possible movements of the bed base; this also applies to stretchers, wheelchairs and chairs where certain parts can be modified or tipped. Check the straps every time the position changes. The abdominal strap attached to the bed frame should compress the mattress.
15. When using the abdominal belt, always use the pelvic belt RFP22000 and lateral belts RFE24000 unless otherwise specified by the prescribing doctor.
16. Ensure that the RENOL key, part number: RFK00100 (RL) and RFK01111 (M) is always available. Use only RENOL locks.
17. In an emergency, if the key is not available or problem with lock system, cut the strap close to the frame.
18. Before use, check operation of the magnet (M) lock system, check correct operation of the RFK01202 Button with the key RFK01111 and pin RFK01201. We recommend replacing the Button RFK01202 after 2 years. The lifetime may be shorter or longer, depending on use.
19. Before use, check operation of the mechanical lock system (RL), check correct operation of the RFK00400 Button with the key RFK00100 and long pin RFK00300. We recommend replacing the Button RFK00400 every 3 years. The lifetime may be shorter or longer, depending on use.
20. The magnetic key (M) must be placed 10 cm minimum from a pacemaker or any other medical system that could be disturbed by a magnet.
21. The magnetic lock system could be activated with a magnet other than that delivered by RENOL.

Maintenance and materials

22. Follow the washing instructions marked on the products.
23. Close Velcro fastenings before washing and drying to prevent them wearing.
24. Detergents: avoid using chlorine.
25. Temperature: maximum temperature 95°C.
26. Drying: delicate cycle at low temperature, about 5% shrinkage possible.
27. Remove all locking systems before washing and drying.
28. Never immerse the magnetic (M) and mechanical closures (RL) in a liquid.
29. Material: 100% Terylene (polyester).
30. The raw materials are treated with fire retardant and anti bacterial.
31. 2 years warranty, this warranty does not cover damage resulting from normal wear and tear.

All our products are CE-approved - Class I medical devices.

Any serious incident occurring in connection with this device must be notified to the manufacturer and to the competent authority of the country where the serious incident took place.

For more information, contact your support specialist, your distributor or www.renolcare.com