



RENOL Positioning and Comfort Solutions

Offer your patients **superior comfort** and **secure positioning** thanks to the combination of the **RCM35209** and **RCV07300** cushions. Developed for care environments, geriatric institutions, and wheelchair users, these aids ensure daily **stability, comfort, and safety**.



RCM35209 – Crescent Cushion

Support, comfort, and stability in a sitting position

The RCM35209 crescent cushion offers exceptional morphological support thanks to its innovative polystyrene micro-granule filling. Designed to perfectly adapt to the shape of the forearms, it ensures optimal ergonomic positioning and effectively prevents slipping.

Its waterproof polyurethane cover guarantees impeccable hygiene, while the adjustable strap with buckle ensures attachment to all types of chairs.



Morphological Comfort

Polystyrene micro-granules that adapt to the shape of the forearms



Ergonomic Positioning

Support for the upper limbs and prevention of slipping



Polyurethane Cover

Waterproof and easy to clean



Adjustable Strap

Attachment to a wheelchair or geriatric chair



- ✓ Ideal for improving the posture and comfort of patients sitting for extended periods
- ✓ Suitable for all types of chairs


Please note: • Once the buckle is closed, if the patient does not have the necessary dexterity to open the strap, the cushion may be considered a restraint. However, this is not its intended use or purpose. For situations requiring restraint, please refer to our Belt range, which is specifically designed for this purpose.



RCV07300 – Visco-elastic Seat Cushion with Abduction Block

Postural Support and Prevention of Sliding

The RCV07300 seat cushion features a **high-density visco-elastic foam** that perfectly conforms to body contours for optimal pressure relief. The **integrated abduction block** effectively prevents forward sliding and promotes a correct and natural posture.



Visco-elastic Foam

Pressure relief and lasting comfort for users who sit for long periods

Integrated Abduction Block

Prevents forward sliding and ensures optimal pelvic positioning; our block is filled with polystyrene beads for added comfort.

Waterproof PU Cover

Maximum hygiene in hospital environments

Universal Compatibility

Fits most wheelchairs and geriatric chairs on the market

✓ Lasting Comfort

Significant reduction in fatigue

✓ Pressure Reduction

Prevention of pressure sores

✓ Stable Posture

Natural and ergonomic support

An ideal combination for global comfort



RCM35209

Optimal support for the upper limbs and forearms



RCV07300

Pelvic stability and postural support



Complete solution

Comfort, safety, and optimal positioning

By combining the **crescent cushion** with the **visco-elastic seat cushion**, you get a **complete positioning solution** that meets the most demanding needs of healthcare professionals and their patients.



Pressure point reduction

The combination of the two cushions ensures homogeneous pressure distribution across the entire body, reducing the risk of discomfort and pressure ulcers.



Prevention of sliding


The abduction block and the adjustable strap work together to keep the patient in a safe and stable position throughout the day.



Posture improvement

Ergonomic positioning of the pelvis and upper limbs promotes a natural posture and reduces muscle fatigue.



 **Ideal for users requiring long-term seated support** – These aids are particularly suitable for long-term care environments, geriatric institutions, and wheelchair users.

Please note: • Once the buckle is closed, if the patient does not have the necessary dexterity to open the strap, the cushion may be considered a restraint. However, this is not its intended use or purpose. For situations requiring restraint, please refer to our Belt range, which is specifically designed for this purpose.

Technical Specifications

Discover the detailed specifications of our RENOL positioning solutions, designed to meet the demands of healthcare professionals.

Reference	Description	Composition	Dimensions	Maintenance
RCM35209	Crescent cushion with belt	PU cover / PS microbeads	130cm	Clean by sponging
RCV07300	Visco seat cushion + abduction cushion	PU cover / Visco-elastic foam	42cmx42cm Seat tilt 3° H- 6cm / 2cm	Clean by sponging
RCV07400	Visco seat cushion + abduction cushion	PU cover / Visco-elastic foam	42cmx42cm Seat tilt 3° H- 10cm / 4cm	Clean by sponging